

# QUICKSTART GUIDE

# PREvention<sup>®</sup> HLD8

## High Level Tool Disinfectant Solution

## Where to Use Prevention<sup>™</sup> HLD8

Prevention<sup>™</sup> HLD8 is a high-level disinfectant for tools and implements. It destroys bacteria, viruses and fungi in 8 minutes!

Perfect for use on:

- Reusable Stainless Steel and Hard Plastic Tools
- Nail Clippers
- Cuticle Scissors and Pushers
- Reusable Drill Bits
- Metal Foot Files
- Comedone Extractors
- Tweezers


### Contact Time

**8 Minutes**  
Bacteria, Viruses & Fungi

Tools and Implements should be fully immersed for the entire contact time.

## How to Use Prevention<sup>™</sup> HLD8


Tool disinfection should take place after each and every client service to prevent the spread of germs. It's important to only use high-quality stainless steel and hard plastic reusable professional tools to avoid compatibility issues.



### STEP ONE

#### Clean, rinse and dry

Clean tools using water and detergent. Rinse and dry thoroughly.



### STEP TWO

#### Fill soaker

Pour Prevention<sup>™</sup> HLD8 directly into soaker. The solution is ready-to-use and should not be diluted.



### STEP THREE

#### Soak tools and re-use solution

Fully immerse tools in Prevention<sup>™</sup> HLD8 for only 8 minutes. Re-use solution for up to 21 days in a soaking tray.\*



### STEP FOUR

#### Rinse, dry and store

Rinse tools with water and dry manually using a clean cloth or disposable towel. Store in a clean container marked "DISINFECTED".

\*The solution must be discarded after 21 days or sooner as dictated by the Prevention<sup>™</sup> HLD8 Solution Test Strips.

## Success & Troubleshooting Tips

### 8 Minutes is all it Takes.

Prevention<sup>™</sup> HLD8 is not intended for use as a holding solution. Tools and implements should be removed after the 8-minute contact time. Excess exposure to the solution does not lead to improved efficacy and may actually lead to wear and tear on your tools over time.

### In Doubt. Throw Out.

Learn to assess the cloudiness of the solution and use your test strips if unsure the product is still effective during the 21 day re-use period. If the solution is visibly contaminated or the test strips show a failing result, throw the solution out and pour fresh.

### Quality is Key.

High quality stainless steel and hard plastic professional tools will have a longer life span after repeated cleaning and disinfection. Low quality tools and carbide bits can face compatibility issues when soaked in disinfectant solutions.

### Always Ready-To-Use.

The solution does not require dilution or mixing with water, making it simple and saving you time. Simply pour Prevention<sup>™</sup> HLD8 directly into the soaking tray and you are ready to disinfect.

For more information visit:  
[LearnAboutRejuvenate.com](http://LearnAboutRejuvenate.com)

**Virox Technologies Inc.**  
2770 Coventry Road  
Oakville, ON L6H 6R1 CANADA  
Toll Free: 1-800-387-7578