

QUICKSTART GUIDE

REJUVENATE® Concentrate

Where to Use Rejuvenate™ Concentrate

Rejuvenate™ Concentrate is an effective and convenient choice for everyday use. It is ideal for facility disinfection, plus it's economical.

Perfect for use on:

- Hard Non-Porous Surfaces and Equipment
- Workstations
- Countertops
- Chairs and Stools
- Circulating* and Non-Circulating Footbaths
- Kitchen and Staff Rooms
- Other Large Surfaces

Dilution	Contact Time
1:16 or 2oz of concentrate to 32oz of water	5 Minutes Bacteria, Viruses & Fungi Surface should be visibly wet for the entire contact time.

How to Use Rejuvenate™ Concentrate


Disinfection should take place after each and every client service or treatment to prevent the spread of germs.



STEP ONE

Dilute the concentrate solution

Dilute Rejuvenate™ Concentrate at a 1:16 dilution ratio (one part Rejuvenate™ with 16 parts water or 2oz of concentrate to 32oz of water).



STEP TWO

Apply diluted solution on surface until completely wet

Remove any visible dirt or debris from surface, or equipment prior to spraying. Diluted solution can be used for up to 90 days in a spray bottle.



STEP THREE

Allow surface to remain visibly wet for 5 minutes

Allowing surfaces to remain wet for 5 minutes kills viruses, bacteria and fungi per label instructions.



AIR DRY

STEP FOUR

Leave to air dry

If desired, you may wipe excess solution using a clean cloth after the required contact time has been achieved.

Success & Troubleshooting Tips

Measure. Measure. Measure.

Using a measuring cup while diluting will help ensure your solution is properly mixed and effective.

More is NOT Better.

Follow the dilution ratios and you will conserve product and prevent unnecessary issues. When diluted properly, Rejuvenate™ Concentrate is a cost-effective solution for your everyday disinfection needs.

Start Fresh.

Bottles and secondary containers should not be topped up with fresh diluted solution. Instead, ensure the containers are emptied, cleaned and dried prior to refilling with new solution to ensure it's effective every time. Remember to also place a workplace label on any new bottles, for easy identification.

For more information visit:
LearnAboutRejuvenate.com

Virox Technologies Inc.
2770 Coventry Road
Oakville, ON L6H 6R1 CANADA
Toll Free: 1-800-387-7578